


AVT-1

Anti Vibration Table

For High Magnification Microscopy

Isolate the microscope from vibrations generated from the external environment



Note:
AVT-1 is designed to be easy to use and almost maintenance free

Technical Specifications

AVT-1

Overall dimensions (W x D x H)	1200 x 800 x 800 mm (47.25" x 31.5" x 31.5")
Weight	70 kg (154 lbs)
Material	Powder painted mild steel Stainless steel and glass
Float dimensions	540 x 340 mm (21.26" x 13.39")



Designed in Denmark



Made in the E.U.

Esco Micro Pte. Ltd.

21 Changi South Street 1 • Singapore 486 777
Tel +65 6542 0833 • Fax +65 6542 6920
csis-medical@escoglobal.com • medical.escoglobal.com

Esco Global Offices: Manama, Bahrain | Beijing, China | Chengdu, China | Guangzhou, China | Shanghai, China
Bangalore, India | Delhi, India | Mumbai, India | Banten, Indonesia | Osaka, Japan | Kuala Lumpur, Malaysia
Melaka, Malaysia | Manila, Philippines | Singapore | Midrand, South Africa | Seoul, South Korea | Salisbury, UK
Philadelphia, PA, USA | Hanoi, Vietnam